

VI ByPass Syringe™ Blue

- Designed for use by nurses
- Includes flush-lock option
- A universal collection device
- Combines traditional products into a simple, intuitive closed system

- Designed to streamline blood collections, reduce false positive blood cultures and make blood collection easier from existing lines.

VI ByPass Plus Syringe

19 gauge side channel permits clinical determination of aspiration volume 0-4 mL

VI Velocity Reduction™

17 gauge cannula helps decrease device-related blood hemolysis

Flush-Lock Option

Line Maintenance

VI Sliding Fluid Seal™


Designed to improve sample quality by removing over 99.9% of contaminants

Prime/Purge Chamber

Helps to capture micro bubbles

Recessed Hub

Standard vacutainer longneck bottle/ tube 1-20 mL



SKU# VI-103-LC

- Description: ByPass Syringe
 - Luer Connection
- Unit/Box: 50
- Unit/Case: 200

Additional Details

- Syringe and Luer Applications
- PIV Insertion
- Blood Cultures and Labs
- Existing Line Lab Draws
- Syringe:
 - 1-20ml sizes
 - Aspiration
 - Flushing
- Vacuum:
 - Tubes
 - Bottles

Evaluated and Endorsed

by **ANTT**[®]

(Aseptic Non Touch Technique)

Reduce contamination and ensure best practice

Learn more by visiting
www.vascularintegrity.com.

