

VI ByPass Syringe™ Red

- Designed for use by Phlebotomists
- A universal collection device
- Combines traditional products into a simple, intuitive closed system
- Designed to streamline blood collections, reduce false positive blood cultures and make blood collection easier from existing lines.

VI Velocity Reduction™

17 gauge cannula helps decrease device-related blood hemolysis

VI ByPass Syringe

19 gauge side channels permit clinical determination of aspiration volume 0-4 mL

VI Sliding Fluid Seal™

Designed to improve sample quality by removing over 99.9% of contaminants

Evaluated and Endorsed



(Aseptic Non Touch Technique)

Reduce contamination and ensure best practice



SKU# VI-102-BC

- Description: ByPass Syringe
 - Blood Culture
- Unit/Box: 50
- Unit/Case: 200

Additional Details

- Direct Vacuum
- Diversion
 - Venipuncture
 - PIV insertion
- Vacuum source
 - Tubes
 - Bottles

Learn more by visiting www.vascularintegrity.com.



See how the VI ByPass Syringe can help you enhance your blood collection processes today.