

## VI-101-WM-A

**VI ByPass Syringe™** without tube adapter.

**Intended Use:** Direct Vacuum Blood collection using Wide Mouth bottles and/or tubes.

### INSTRUCTIONS

Prep supplies and site per hospital policy according to Aseptic Non-Touch Technique (ANTT).

**STEP 1:** Attach the VI ByPass Syringe™ to a Butterfly extension or IV extension set. *\*Disinfect bottle tops according to hospital policy.*

**STEP 2:** Cannulate vessel, pull back on the plunger, back prime/aspirate and sequester up to 4mls into the VI ByPass Syringe™.

**STEP 3:** Hold the back portion of the ByPass Syringe and insert/push bottle upward over the booted needle until blood flows into the bottle. Maintain bottles in upright position and allow to fill to the correct volume.

**STEP 4:** Sequestered Blood may be reinfused according to procedural policy by pushing forward on the syringe plunger or discarded within the device.

**STEP 5:** Remove VI ByPass Syringe™ and dispose of sharps container per hospital policy.

Perform remaining Phlebotomy or line maintenance steps per hospital policy.

 REF	- CATALOG NUMBER	 STERILE	 R	- STERILIZED USING IRRADIATION
 LOT	- LOT NUMBER			- DO NOT USE IF PACKAGE IS DAMAGED
	- NON-PYROGENIC			- DO NOT RESTERILIZE
	- CONSULT OPERATING INSTRUCTIONS		Rx only	- RX ONLY
	- DO NOT REUSE			USE BY: xx-xx-xxxx
				- LATEX FREE
	Vascular Integrity 16661 Jamboree Rd. Irvine, CA 92606		ICU Medical 951 Calle Amanecer San Clemente, CA 92673	

## VI-101-WM-A

**VI ByPass Syringe™** without tube adapter.

**Intended Use:** Direct Vacuum Blood collection using Wide Mouth bottles and/or tubes.

### INSTRUCTIONS

Prep supplies and site per hospital policy according to Aseptic Non-Touch Technique (ANTT).

**STEP 1:** Attach the VI ByPass Syringe™ to a Butterfly extension or IV extension set. *\*Disinfect bottle tops according to hospital policy.*

**STEP 2:** Cannulate vessel, pull back on the plunger, back prime/aspirate and sequester up to 4mls into the VI ByPass Syringe™.

**STEP 3:** Hold the back portion of the ByPass Syringe and insert/push bottle upward over the booted needle until blood flows into the bottle. Maintain bottles in upright position and allow to fill to the correct volume.

**STEP 4:** Sequestered Blood may be reinfused according to procedural policy by pushing forward on the syringe plunger or discarded within the device.

**STEP 5:** Remove VI ByPass Syringe™ and dispose of sharps container per hospital policy.

Perform remaining Phlebotomy or line maintenance steps per hospital policy.

 REF	- CATALOG NUMBER	 STERILE	 R	- STERILIZED USING IRRADIATION
 LOT	- LOT NUMBER			- DO NOT USE IF PACKAGE IS DAMAGED
	- NON-PYROGENIC			- DO NOT RESTERILIZE
	- CONSULT OPERATING INSTRUCTIONS		Rx only	- RX ONLY
	- DO NOT REUSE			USE BY: xx-xx-xxxx
				- LATEX FREE
	Vascular Integrity 16661 Jamboree Rd. Irvine, CA 92606		ICU Medical 951 Calle Amanecer San Clemente, CA 92673	