

VI ByPass Syringe™

Blue Device with back-end luer access
VI-103-LC

Intended Use:

VI ByPass Syringe™ LC is compatible with Syringes and vacuum holders for blood collection and existing line procedures that require flushing.

Directions for Use

STEP 1: Prep supplies- Following **Standard Aseptic Non-Touch Technique** attach the VI ByPass Syringe™ (Blue) to a Butterfly extension, IV extension set or central line with or without needleless connector. An extension set must be used when using new or existing peripheral IV catheter.

STEP 2: Aspirate and sequester up to 4mls per policy for Blood or Blood Culture collection. The ByPass syringe will prime during aspiration. Collect Blood via Syringe or attach vacuum holder for sample collection.

STEP 3: Remove vacuum holder device from back-end luer.

STEP 4: Reinfuse sequestered volume or not per facility policy.

STEP 5: Attach prefilled saline syringe, Flush and Lock per facility policy

STEP 6: Detach flush syringe and VI ByPass Syringe™ from the line following facility policy for clamping and discard per hospital protocol.

 REF	- CATALOG NUMBER	 STERILE	R	- STERILIZED USING IRRADIATION
 LOT	- LOT NUMBER			- DO NOT USE IF PACKAGE IS DAMAGED
	- NON-PYROGENIC			- DO NOT RESTERILIZE
	- CONSULT OPERATING INSTRUCTIONS		Rx only	- RX ONLY
	- DO NOT REUSE			USE BY: xx-xx-xxxx
	Vascular Integrity 16661 Jamboree Rd. Irvine, CA 92606		ICU Medical 951 Calle Amanecer San Clemente, CA 92673	- LATEX FREE

VI ByPass Syringe™

Blue Device with back-end luer access
VI-103-LC

Intended Use:

VI ByPass Syringe™ LC is compatible with Syringes and vacuum holders for blood collection and existing line procedures that require flushing.

Directions for Use

STEP 1: Prep supplies- Following **Standard Aseptic Non-Touch Technique** attach the VI ByPass Syringe™ (Blue) to a Butterfly extension, IV extension set or central line with or without needleless connector. An extension set must be used when using new or existing peripheral IV catheter.

STEP 2: Aspirate and sequester up to 4mls per policy for Blood or Blood Culture collection. The ByPass syringe will prime during aspiration. Collect Blood via Syringe or attach vacuum holder for sample collection.

STEP 3: Remove vacuum holder device from back-end luer.

STEP 4: Reinfuse sequestered volume or not per facility policy.

STEP 5: Attach prefilled saline syringe, Flush and Lock per facility policy

STEP 6: Detach flush syringe and VI ByPass Syringe™ from the line following facility policy for clamping and discard per hospital protocol.

 REF	- CATALOG NUMBER	 STERILE	R	- STERILIZED USING IRRADIATION
 LOT	- LOT NUMBER			- DO NOT USE IF PACKAGE IS DAMAGED
	- NON-PYROGENIC			- DO NOT RESTERILIZE
	- CONSULT OPERATING INSTRUCTIONS		Rx only	- RX ONLY
	- DO NOT REUSE			USE BY: xx-xx-xxxx
	Vascular Integrity 16661 Jamboree Rd. Irvine, CA 92606		ICU Medical 951 Calle Amanecer San Clemente, CA 92673	- LATEX FREE